

Product Datasheet

FGF21 Antibody / Fibroblast Growth Factor 21 (orb699561)

Description

Fibroblast growth factor-1 (FGF-1), also designated


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P

IHC staining of FFPE human kidney with F...

Immunogen

A portion of amino acids 91-209 from the human protein was used as the immunogen for the FGF21 antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the FGF21 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the FGF21 antibody should be determined by the researcher.

IHC staining of FFPE human breast with F...

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2b, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID
[Q9NSA1](#)
Hazard Information

This FGF21 antibody is available for research use only.


Dilution Range

Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT

Expiration Date

12 months from date of receipt.

Human Protein Microarray Specificity Validation



Analysis of HuProt (TM) microarray conta...