
Product Datasheet

**TTR Antibody / Transthyretin / Prealbumin
(orb699586)**

Description

Prealbumin, also designated transthyretin, is a major

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

IHC-P, WB

Immunogen

Recombinant human full-length protein was used as the immunogen for the Transthyretin antibody.

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the Transthyretin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the Transthyretin antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2c, kappa

Clonality

Monoclonal

Purity

Protein G affinity chromatography

Uniprot ID

P02766

Hazard Information

This Transthyretin antibody is available for research use only.

Dilution Range


Western blot: 2-4ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT

Expiration Date


12 months from date of receipt.



Western blot testing of human liver lysa...



IHC staining of FFPE human liver with Tr...



Analysis of HuProt (TM) microarray conta...