



Product Datasheet

Factor VII Antibody (orb699636)


Description	Hemostasis following tissue injury involves the deployment of essential plasma procoagulants (prothr...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	A portion of amino acids 366-466 from the human protein was used as the immunogen for the Factor VII antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the Factor VII antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the Factor VII antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P08709
Hazard Information	This Factor VII antibody is available for research use only.
Dilution Range	ELISA: order Ab without BSA for coating, Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.



IHC staining of FFPE human liver with Fa...



IHC staining of FFPE human liver with Fa...



Analysis of HuProt (TM) microarray conta...