


Product Datasheet


HSP90 beta Antibody / HSP90AB1 (orb699673)

Description	This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signa...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	A portion of amino acids 581-704 from the human protein was used as the immunogen for the HSP90AB1 antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the HSP90AB1 antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the HSP90AB1 antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG1, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P08238
Hazard Information	This HSP90AB1 antibody is available for research use only.
Dilution Range	Western blot: 1-2ug/ml, Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.


Human Protein Microarray Specificity Validation



Analysis of HuProt (TM) microarray conta...



IHC staining of FFPE human pancreas with...



Western blot testing of human MCF7 cell ...