


Product Datasheet


IL1 alpha Antibody / IL1A (orb699766)

Description	Two forms of interleukin-1, designated IL-1 alpha and beta, have been described. Although encoded by...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	A portion of amino acids 113-271 from the human protein was used as the immunogen for the IL1A antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Store the IL1A antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).
Note	For research use only
Application notes	Optimal dilution of the IL1A antibody should be determined by the researcher.
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Purity	Protein G affinity chromatography
Uniprot ID	P01583
Hazard Information	This IL1A antibody is available for research use only.
Dilution Range	Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT
Expiration Date	12 months from date of receipt.


Human Protein Microarray Specificity Validation



Analysis of HuProt (TM) microarray conta...



IHC staining of FFPE human colon carcino...



SDS-PAGE analysis of purified, BSA-free ...