
Product Datasheet

TG Antibody / Thyroglobulin (orb699783)

Description

mAb rTGB/4744 reacts with a partially defined epitope


Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested Applications

FACS, IHC-P

Immunogen

Human thyroid follicular cells were used as the immunogen for the recombinant Thyroglobulin antibody.

IHC staining of FFPE human thyroid carci...

Preservatives

0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide

Storage

Store the recombinant Thyroglobulin antibody at 2-8°C (with azide) or aliquot and store at -20°C or colder (without azide).

Note

For research use only

Application notes

Optimal dilution of the recombinant Thyroglobulin antibody should be determined by the researcher.

Formula

0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide

Isotype

Mouse IgG2a, kappa

Clonality

Recombinant

Purity

Protein G affinity chromatography

Uniprot ID

P01266

Hazard Information

This recombinant Thyroglobulin antibody is available for research use only.

Dilution Range

Flow cytometry: 0.5-1ug/million cells, Immunohistochemistry (FFPE): 1-2ug/ml for 30 minutes at RT

Expiration Date

12 months from date of receipt.