

Product Datasheet

Goat anti-Rabbit IgG (H&L) antibody (Biotin) (orb699949)

Description	Goat polyclonal antibody against Goat anti-Rabbit IgG (H&L) conjugated to Biotin
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Purified Rabbit IgG, whole molecule
Form/Appearance	Liquid
Storage	Store non-diluted antibody at 2-8°C. Shelf life is 1 year from date of receipt.
Note	For research use only
Application notes	No reactivity is observed to non-immunoglobulin rabbit serum proteins based in immunoelectrophoresis
Clonality	Polyclonal
Purity	Immunogen purified goat IgG.
Dilution Range	1 : 2000-1:20 000 (ELISA), 1:200 (ICC), 1 : 500 -1 : 5000 (IHC), 1 : 2000-1:20 000 (WB)
Expiration Date	12 months from date of receipt.