
Product Datasheet

**Donkey anti-Sheep IgG (H&L) antibody (ALP)
(orb699955)**

Description

Donkey polyclonal antibody against

conjugated to AP

Species/Host

Donkey

Reactivity

Sheep

Conjugation

ALP

Tested Applications

ELISA, ICC, IHC, WB

Immunogen

Purified sheep IgG, whole molecule

Form/Appearance

Liquid

Storage

Non-diluted antibody is stable for 4 years at 2-8°C. For storage at -20°C dilute antibody solution with an equal volume of glycerol to obtain final glycerol concentration of 50 % to prevent loss of enzymatic activity. Such solution will not freeze in -20°C. If you are using a 1:5000 dilution prior to diluting with glycerol, then you would need to use a 1:2500 dilution after adding glycerol. Prepare working dilution prior to use and then discard, Be sure to mix well but without foaming.

Note

For research use only

Application notes

No reactivity is observed to non-immunoglobulin sheep serum proteins based in immunoelectrophoresis. BSA and milk have to be replaced by other blocking reagents, like donkey serum or commercial formulations which are free from bovine IgG.

Clonality

Polyclonal

Purity

Immunogen affinity purified donkey IgG.

Dilution Range

1 : 50-1 : 5000 (ICC), 1 : 500-1 : 2000 (ELISA), 1 : 20 -1 : 2000 (IHC), 1 : 500-1 : 2000 (WB)

Expiration Date

12 months from date of receipt.