

Product Datasheet

Goat anti-Rabbit IgG (H&L) antibody (FITC) (orb699967)

Description	Goat polyclonal antibody against Goat anti-Rabbit IgG (H&L) conjugated to FITC
Species/Host	Goat
Conjugation	FITC
Form/Appearance	Liquid
Storage	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
Note	For research use only
Application notes	This antibody reacts with the heavy chains on rabbit IgG and with the light chains on all rabbit immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin rabbit serum proteins based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1% (w/v) BSA, protease/IgG free and 0.05 % (w/v) sodium azide as preservative.
Clonality	Polyclonal
Purity	Antigen Immunogen affinity purified IgG, F(ab)' ₂ fragment.
Dilution Range	1 : 50-1 : 5000 (ICC) (Flow cyt), 1 : 20-1 : 2000 (IHC)
Expiration Date	12 months from date of receipt.