

## Product Datasheet

# Donkey anti-Sheep IgG (H&L) antibody (FITC) (orb699979)

<b>Description</b>	Donkey polyclonal antibody against Donkey anti-Sheep IgG (H&L) conjugated to FITC
<b>Species/Host</b>	Donkey
<b>Conjugation</b>	FITC
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	This antibody reacts with the heavy chains on sheep IgG and with the light chains on all sheep immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin sheep serum proteins and human, mouse, rabbit IgG based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7, 1% (w/v) BSA, Protease/IgG free and 0.05 % (w/v) sodium azide as preservative. BSA and milk have to be replaced by other blocking reagents, like donkey serum or commercial formulations which are free from bovine IgG.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Immunogen affinity purified IgG.
<b>Dilution Range</b>	1 : 20-1 : 2000 (IHC), 1 : 50-1 : 5000 (ICC) (Flow cyt)
<b>Expiration Date</b>	12 months from date of receipt.