

Product Datasheet

Goat anti-Rat IgG (H&L) antibody (FITC) (orb699990)

| | |
|--------------------------|--|
| Description | Goat polyclonal antibody against Goat anti-Rat IgG (H&L) conjugated to FITC |
| Species/Host | Goat |
| Conjugation | FITC |
| Form/Appearance | Liquid |
| Storage | Store at 2-8°C. Shelf life of this product is one year from date of receipt. |
| Note | For research use only |
| Application notes | This antibody reacts with the heavy chains on rat IgG and with the light chains on all rat immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin rat serum proteins or to human or mouse IgG based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1% (w/v) BSA, protease/IgG free and 0.05 % (w/v) sodium azide as preservative. |
| Clonality | Polyclonal |
| Purity | Immunogen affinity purified IgG. |
| Dilution Range | 1 : 50-1 : 5000 (ICC), 1 : 20-1 : 2000 (IHC) |
| Expiration Date | 12 months from date of receipt. |