

Product Datasheet

Goat anti-Rabbit IgG (H&L antibody (Biotin) (orb699999)

Description	Goat polyclonal antibody against Goat anti-Rabbit IgG (H&L conjugated to Biotin)
Species/Host	Goat
Conjugation	Biotin
Form/Appearance	Liquid
Storage	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
Note	For research use only
Application notes	<p>This antibody reacts with the heavy chains on rabbit IgG and with the light chains on all rabbit immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin rabbit serum proteins based on immunoelectrophoresis. Minimum cross-reactivity is observed to IgG from bovine, goat, human, mouse or rat based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free and 0.1 % (v/v) Kathon CG is used as preservative. Use of sodium azide will inhibit enzymatic activity of horseradish peroxidase.</p>
Clonality	Polyclonal
Purity	Immunogen affinity purified IgG.
Dilution Range	The optimal working dilution should be determined by the investigator
Expiration Date	12 months from date of receipt.