

Product Datasheet

Goat anti-Rat IgG (H&L) antibody (TRITC) (orb700021)

Description	Goat polyclonal antibody against Goat anti-Rat IgG (H&L) conjugated to TRITC
Species/Host	Goat
Conjugation	TRITC
Form/Appearance	Liquid
Storage	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
Note	For research use only
Application notes	This antibody reacts with the heavy chains on rat IgG and with the light chains on all rat immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin rat serum proteins or IgG from human or mouse based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1% (w/v) BSA, protease/IgG free and 0.05 % (w/v) sodium azide as preservative.
Clonality	Polyclonal
Purity	Immunogen affinity purified IgG, F(ab)' ₂ fragment.
Dilution Range	The optimal working dilution should be determined by the investigator, (ICC), immunohistochemistry (IHC) (Flow cyt)
Expiration Date	12 months from date of receipt.