

## Product Datasheet

# Donkey anti-Rat IgG (H&L antibody (Biotin) (orb700026)

<b>Description</b>	Donkey polyclonal antibody against Donkey anti-Rat IgG (H&L conjugated to Biotin
<b>Species/Host</b>	Donkey
<b>Conjugation</b>	Biotin
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	This antibody reacts with the heavy chains on rat IgG and with the light chains on all rat immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin rabbit serum proteins based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free and 0.1 % (v/v) Kathon CG is used as preservative.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Immunogen affinity purified IgG.
<b>Dilution Range</b>	The optimal working dilution should be determined by the investigator
<b>Expiration Date</b>	12 months from date of receipt.