

Product Datasheet

Donkey anti-Rat IgG (H&L) antibody (TRITC) (orb700034)

Description	Donkey polyclonal antibody against Donkey anti-Rat IgG (H&L) conjugated to TRITC
Species/Host	Donkey
Conjugation	TRITC
Form/Appearance	Liquid
Storage	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
Note	For research use only
Application notes	This antibody reacts with the heavy chains on rat IgG and with the light chains on all rat immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin rat serum proteins based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7, 0.05 % (w/v) sodium azide as preservative. 1% (w/v) BSA, Protease/IgG free.
Clonality	Polyclonal
Purity	Immunogen affinity purified IgG.
Dilution Range	The optimal working dilution should be determined by the investigator,
Expiration Date	12 months from date of receipt.