

Product Datasheet

Goat anti-Mouse IgG (H&L) antibody (TRITC) (orb700047)

Description	Goat polyclonal antibody against Goat anti-Mouse IgG (H&L) conjugated to TRITC
Species/Host	Goat
Conjugation	TRITC
Form/Appearance	Liquid
Storage	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
Note	For research use only
Application notes	This antibody reacts with the heavy chains on mouse IgG and with the light chains on all mouse immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin mouse serum proteins and goat, human, rabbit, rat IgG based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1% (w/v) BSA, protease/IgG free and 0.05 % (w/v) sodium azide as preservative.
Clonality	Polyclonal
Purity	Immunogen affinity purified goat IgG.
Dilution Range	The optimal working dilution should be determined by the investigator
Expiration Date	12 months from date of receipt.