

## Product Datasheet

# Goat anti-Rat IgG (H&L) antibody (FITC) (orb700059)

<b>Description</b>	Goat polyclonal antibody against Goat anti-Rat IgG (H&L) conjugated to FITC
<b>Species/Host</b>	Goat
<b>Conjugation</b>	FITC
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	This antibody reacts with the heavy chains on rat IgG and with the light chains on all rat immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin rat serum proteins based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1% (w/v) BSA, protease/IgG free and 0.05 % (w/v) sodium azide as preservative.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Immunogen affinity purified IgG.
<b>Dilution Range</b>	1 : 50-1 : 5000 (ICC) (Flow cyt), 1 : 20-1 : 2000 (IHC)
<b>Expiration Date</b>	12 months from date of receipt.