

---

## Product Datasheet

### Donkey anti-Sheep IgG (H&L antibody (Biotin) (orb700070)

**Description**

Donkey polyclonal antibody against

to Biotin

**Species/Host**

Donkey

**Conjugation**

Biotin

**Form/Appearance**

Liquid

**Storage**

Store at 2-8°C. Shelf life of this product is one year from date of receipt.

**Note**

For research use only

**Application notes**

This antibody reacts with the heavy chains on sheep IgG and with the light chains on all sheep immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin sheep serum proteins and human and rabbit IgG based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free and 0.1 % (v/v) Kathon CG is used as preservative. BSA and milk have to be replaced by other blocking reagents, like donkey serum or commercial formulations which are free from bovine IgG.

**Clonality**

Polyclonal

**Purity**

Immunogen affinity purified IgG.

**Dilution Range**

The optimal working dilution should be determined by the investigator

**Expiration Date**

12 months from date of receipt.