

Product Datasheet

Donkey anti-Sheep IgG (H&L) antibody (FITC) (orb700071)

Description	Donkey polyclonal antibody against Donkey anti-Sheep IgG (H&L) conjugated to FITC
Species/Host	Donkey
Conjugation	FITC
Form/Appearance	Liquid
Storage	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
Note	For research use only
Application notes	This antibody reacts with the heavy chains on sheep IgG and with the light chains on all sheep immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin sheep serum proteins and human, rabbit IgG based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7, 1% (w/v) BSA, Protease/IgG free, 0.05 % (w/v) sodium azide as preservative. BSA and milk have to be replaced by other blocking reagents, like donkey serum or commercial formulations which are free from bovine IgG.
Clonality	Polyclonal
Purity	Immunogen affinity purified IgG.
Dilution Range	1 : 20-1 : 2000 (IHC), 1 : 50-1 : 5000 (ICC) (Flow cyt)
Expiration Date	12 months from date of receipt.