
Product Datasheet

**Donkey anti-Sheep IgG (H&L) antibody (FITC)
(orb700087)**

Description

Donkey polyclonal antibody against

to FITC

Species/Host

Donkey

Conjugation

FITC

Form/Appearance

Liquid

Storage

Store at 2-8°C. Shelf life of this product is one year from date of receipt.

Note

For research use only

Application notes

This antibody reacts with the heavy chains on sheep IgG and with the light chains on all sheep immunoglobulins based on immunoelectrophoresis. Minimum cross-reactivity is observed to non-immunoglobulin sheep serum proteins based on immunoelectrophoresis. Antibody is supplied in 10 mM sodium phosphate, 150 mM sodium chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free and 0.1 % (v/v) Kathon CG is used as preservative. BSA and milk have to be replaced by other blocking reagents, like donkey serum or commercial formulations which are free from bovine IgG.

Clonality

Polyclonal

Purity

Immunogen affinity purified IgG.

Dilution Range

1 : 20-1 : 2000 (IHC), 1 : 50-1 : 5000 (ICC) (Flow cyt)

Expiration Date

12 months from date of receipt.