

## Product Datasheet

### Donkey anti-Goat IgG (H&L) antibody (Biotin) (orb700175)

<b>Description</b>	Donkey polyclonal antibody against Donkey anti-Goat IgG (H&L) conjugated to Biotin
<b>Species/Host</b>	Donkey
<b>Conjugation</b>	Biotin
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	This antibody reacts with the heavy chains on goat IgG and with the light chains on all goat immunoglobulins based on immunoelectrophoresis.No reactivity is observed to non-immunoglobulin goat serum proteins based on immunoelectrophoresis.BSA and milk have to be replaced by other blocking reagents, like donkey serum or commercial formulations which are free from bovine IgG.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Immunogen affinity purified IgG.
<b>Dilution Range</b>	The optimal working dilution should be determined by the investigator
<b>Expiration Date</b>	12 months from date of receipt.