

Product Datasheet

Goat anti-Human IgA + IgG + IgM antibody (TRITC) (orb700176)

Description	Goat polyclonal antibody against Goat anti-Human IgA + IgG + IgM conjugated to TRITC
Species/Host	Goat
Conjugation	TRITC
Form/Appearance	Liquid
Storage	Store at 2-8°C. Shelf life of this product is one year from date of receipt.
Note	For research use only
Application notes	This antibody reacts with the heavy and light chains on human IgG, IgA and IgM based on immunoelectrophoresis. No reactivity is observed to non-immunoglobulin human serum proteins based on immunoelectrophoresis.
Clonality	Polyclonal
Purity	Immunogen affinity purified IgG, F(ab)' ₂ fragment.
Dilution Range	The optimal working dilution should be determined by the investigator,
Expiration Date	12 months from date of receipt.