

Product Datasheet

Goat anti-Guinea pig IgG (H&L) antibody (Biotin) (orb700278)

Description	Goat polyclonal antibody against Goat anti-Guinea pig IgG (H&L) conjugated to Biotin
Species/Host	Goat
Conjugation	Biotin
Tested Applications	ELISA, IHC, WB
Immunogen	Purified Guinea pig IgG, whole molecule
Form/Appearance	Liquid
Storage	Store non-diluted antibody at 2-8°C. Shelf life is 1 year from date of receipt.
Note	For research use only
Application notes	No reactivity is observed to non-immunoglobulin guinea pig serum proteins based in immunoelectrophoresis. This antibody reacts with the heavy chains on guinea pig IgG and with the light chains on all guinea pig immunoglobulins based on immunoelectrophoresis.
Clonality	Polyclonal
Purity	Immunogen affinity purified goat IgG.
Dilution Range	WB: 1: 2000-20000, IHC-P: 1: 500-1: 5000
Expiration Date	12 months from date of receipt.