

Product Datasheet

Rabbit anti-Goat IgG (H&L), Biotin conjugated, min, cross-reactivity to human, mouse and Rat IgG (orb700281)

Description	Rabbit polyclonal antibody against Rabbit anti-Goat IgG (H&L conjugated to Biotin
Species/Host	Rabbit
Reactivity	Goat
Conjugation	Biotin
Tested Applications	ELISA, IHC, WB
Immunogen	Purified Goat IgG, whole molecule
Form/Appearance	Liquid
Storage	Store non-diluted antibody at 2-8°C. Shelf life is 1 year from date of receipt.
Note	For research use only
Application notes	No reactivity is observed to non-immunoglobulin goat serum proteins based in immunoelectrophoresis.No reactivity is observed to human, mouse and rat IgG based on immunoelectrophoresis.
Clonality	Polyclonal
Purity	Immunogen affinity purified rabbit IgG.
Dilution Range	1 : 2000-20 000 (ELISA), 1 : 500 -1 : 5000 (IHC), 1 : 2000-20 000 (WB)
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)