

Product Datasheet

Goat anti-Rabbit IgG Fc antibody (Biotin) (orb700368)

Description	Goat polyclonal antibody against Goat anti-Rabbit IgG Fc conjugated to Biotin
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	IA
Immunogen	Purified Rabbit IgG
Form/Appearance	Liquid
Storage	Store non-diluted antibody at 2-8°C. Shelf life is 1 year from date of receipt.
Note	For research use only
Application notes	This antibody reacts with the heavy chains on rabbit IgG based on immunoelectrophoresis. No reactivity is observed to immunoglobulin light chains or non-immunoglobulin rabbit serum proteins based on immunoelectrophoresis.
Clonality	Polyclonal
Purity	Immunogen affinity purified goat IgG.
Dilution Range	The optimal working dilution should be determined by the investigator
Expiration Date	12 months from date of receipt.