
Product Datasheet

**Rabbit anti-Bovine IgG (H&L) antibody (HRP)
(orb700453)**

Description

Rabbit polyclonal antibody against

conjugated to HRP

Species/Host

Rabbit

Reactivity

Bovine

Conjugation

HRP

Tested Applications

ELISA, IHC, WB

Immunogen

Purified bovine IgG, whole molecule

Form/Appearance

Lyophilized

Storage

Store freeze-dried powder at 2-8 °C. For long term storage after reconstitution, dilute with 50% glycerol and store at -20 °C as a liquid.

Note

For research use only

Application notes

Antibody purity is > 95% based on SDS-PAGE. Antibody is supplied in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free. 0.1% (v/v) ProClin 150 is added as a preservative.

Clonality

Polyclonal

Purity

Immunogen affinity purified rabbit IgG.

Uniprot ID

DAA31645

Dilution Range

1 : 200-1 : 5 000 (ELISA), 1 : 500-1 : 5 000 (IHC), 1 : 200-1 : 5 000 (WB)

Expiration Date

12 months from date of receipt.