
Product Datasheet

**Goat anti-Mouse IgG (H&L) antibody (TRITC)
(orb700537)**

Description

Goat polyclonal antibody against Goat

	TRITC
Species/Host	Goat
Reactivity	Mouse
Conjugation	TRITC
Tested Applications	FC, IL
Immunogen	Purified Mouse IgG, whole molecule
Form/Appearance	Lyophilized
Storage	Store at -20° C.
Note	For research use only
Application notes	Antibody purity is > 90% based on SDS-PAGE Small amounts of intact IgG may be present.Affinity purified using solid phase Mouse IgG .Antibody is supplied in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free.Based on IEP, this antibody reacts with: heavy (γ) chains on mouse IgG light chains on all mouse immunoglobulins .Based on IEP, no reactivity is observed to: non-immunoglobulin human serum immunoglobulins bovine, goat, human, rabbit or rat IgG .
Clonality	Polyclonal
Purity	Immunogen affinity purified goat IgG.
Dilution Range	1 : 20-1 : 2000 (Flow cyt), (IL)
Expiration Date	12 months from date of receipt.