
Product Datasheet

**Goat anti-Mouse IgG + IgM antibody (HRP)
(orb700542)**

Description

Goat polyclonal antibody against Goat

	HRP
Species/Host	Goat
Reactivity	Mouse
Conjugation	HRP
Tested Applications	FC, IL
Immunogen	Purified Mouse IgG + IgM, whole molecule
Form/Appearance	Lyophilized
Storage	Store freeze-dried powder at 2-8 °C.
Note	For research use only
Application notes	Antibody purity is > 95% based on SDS-PAGE. Affinity purified using solid phase Mouse IgG + IgM. Antibody is supplied in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free. 0.1% (v/v) Kathon CG is added as a preservative. Based on IEP, this antibody reacts with: heavy chains on mouse IgG + IgM light chains on all mouse immunoglobulins. Based on IEP, no reactivity is observed to: non-immunoglobulin mouse serum proteins.
Clonality	Polyclonal
Purity	Immunogen affinity purified goat IgG.
Dilution Range	1 : 20-1 : 2000 (Flow cyt), (IL)
Expiration Date	12 months from date of receipt.