

---

## Product Datasheet

**Goat anti-Mouse IgM (Å,Âµ chain) antibody (HRP)  
(orb700555)**

**Description**

Goat polyclonal antibody against

conjugated to HRP

**Species/Host**

Goat

**Reactivity**

Mouse

**Conjugation**

HRP

**Tested Applications**

FC, IL

**Immunogen**

 Purified Mouse IgM,  $\mu$  chain

**Form/Appearance**

Lyophilized

**Storage**

Store freeze-dried powder at 2-8 °C.

**Note**

For research use only

**Application notes**

Antibody purity is > 95% based on SDS-PAGE. Affinity purified using solid phase Mouse IgM. Antibody is supplied in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 1 % (w/v) BSA, Protease/IgG free. 0.1% (v/v) Kathon CG is added as a preservative. Based on IEP, this antibody reacts with: heavy ( $\gamma$ ) chains on mouse IgM. Based on IEP, no reactivity is observed to: non-immunoglobulin mouse serum proteins light chains on all mouse immunoglobulins.

**Clonality**

Polyclonal

**Purity**

Immunogen affinity purified goat IgG.

**Uniprot ID**

AAA64614

**Dilution Range**

1 : 20-1 : 2000 (Flow cyt), (IL)

**Expiration Date**

12 months from date of receipt.