

## Product Datasheet

### Human Syndecan-1 Protein (orb700634)

<b>Description</b>	Human Syndecan-1, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 23 - Gly 254...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Syndecan-1
<b>Preservatives</b>	PBS, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-Fc
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 50.3 kDa. The protein migrates as 65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	Gln 23 - Gly 254
<b>Purity</b>	95%
<b>Source</b>	Human Syndecan-1, Fc Tag (orb700634) is expressed from human 293 cells (HEK293). It contains AA Gln 23 - Gly 254 (Accession # P18827-1).
<b>MW</b>	50.3 kDa
<b>Uniprot ID</b>	<a href="https://www.uniprot.org/entry/P18827">P18827</a>
<b>Expiration Date</b>	6 months from date of receipt.