
Product Datasheet

Biotinylated Human Syndecan-1 Protein (orb700635)

Description

Biotinylated Human Syndecan-1, His, is

(HEK293). It contains AA Gln 2...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

Syndecan-1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His & C-Avi

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (.). The protein has a calculated MW of 27.5 kDa. The protein migrates as 44-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

Gln 23 - Gly 254

Purity

90%

Source

Biotinylated Human Syndecan-1, His, (orb700635) is expressed from human 293 cells (HEK293). It contains AA Gln 23 - Gly 254 (Accession # P18827-1).

MW

27.5 kDa

Uniprot ID**P18827****Expiration Date**

6 months from date of receipt.