

Product Datasheet

Cynomolgus Siglec-10 Protein (orb700637)

Description	Cynomolgus Siglec-10, His Tag is expressed from human 293 cells (HEK293). It contains AA Thr 17 - Al...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Siglec-10
Preservatives	25 mM MES, 150 mM NaCl, 0.5 M Arginine, pH5.0
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 60.2 kDa. The protein migrates as 65-75 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Thr 17 - Ala 549
Purity	95%
Source	Cynomolgus Siglec-10, His Tag (orb700637) is expressed from human 293 cells (HEK293). It contains AA Thr 17 - Ala 549 (Accession # A0A2K5WBX8-1).
MW	60.2 kDa
Uniprot ID	A0A2K5WBX8
Expiration Date	6 months from date of receipt.