



Product Datasheet

CD9 antibody (orb704168)


Description	CD9 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human CD9 (111-196/228 aa)
Target	CD9
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	25 kDa
Uniprot ID	P21926
Dilution Range	WB=1:500-1000, IHC-P=1:50-200, IHC-F=1:50-200, ICC=1:100, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P analysis of human laryngeal carcin...



IHC-P analysis of mouse kidney using CD9...



Western blot: Lane 1:Kidney (Mouse);Lane

...