



Product Datasheet

Gelsolin antibody (orb704177)


Description	Gelsolin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human Gelsolin protein
Target	Gelsolin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	80 kDa
Uniprot ID	P06396
Dilution Range	WB=1:500-2000, IHC-P=1:50-200
Expiration Date	12 months from date of receipt.



IHC-P analysis of mouse pancreas using G...



Western blot: Lane 1: Plasma (Mouse); Lane ...



Western blot: Lane 1: HepG2 (Human) Cell; L...