

Product Datasheet

Syk (phospho-Tyr323) antibody (orb704290)

Description	Syk (phospho-Tyr323) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Syk around the phosphorylation site of Tyr323 (NP(p-Y)EP)
Target	Phospho-Syk (Tyr323)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	72 kDa
Uniprot ID	P43405
Expiration Date	12 months from date of receipt.