



Product Datasheet

NADPH Oxidase 4 antibody (orb704314)


Description	NADPH Oxidase 4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	Recombinant human NADPH oxidase 4 (C-terminal)
Target	NADPH oxidase 4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	67 kDa
Uniprot ID	Q9NPH5
Dilution Range	WB=1:1000, IHC-P=1:50-200, ICC=1:50
Expiration Date	12 months from date of receipt.



IHC-P analysis of human lung using NADPH...



IHC-P analysis of kidney using NADPH oxi...



IHC-P analysis of mouse lung using NADPH...