


Product Datasheet

Nano-Tag (15) antibody (orb704437)

| | |
|----------------------------|---|
| Description | Nano-Tag (15) antibody |
| Species/Host | Rabbit |
| Reactivity | All |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | KLH conjugated synthetic peptide Nano-Tag15 |
| Target | Nano-Tag (15) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/1ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Dilution Range | WB=1:1000-2000 |
| Expiration Date | 12 months from date of receipt. |



Western blot: Lane 1: Nano-Tagged Fusion P...