


Product Datasheet

BCMA Antibody (orb704540)

Description	Anti-BCMA (belantamab biosimilar) mAb
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	BCMA
Form/Appearance	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	IgG1
Clonality	Recombinant
Purity	Purified from cell culture supernatant by affinity chromatography
Uniprot ID	Q02223
Dilution Range	ELISA 1:5000-10000; Flow Cyt 1:101
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...