

Product Datasheet

Human Transmembrane protease serine 2 protein (orb704556)

Description	Recombinant Human Transmembrane protease serine 2(TMPRS2),partial
Endotoxins	Not test.
Preservatives	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal GST-tagged
Note	For research use only
Application notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Protein Sequence	IVGGESALPGAWPWQVSLHVQNVHVCGGSITPEWIVTAAHCVEKPLNPPWHW TAFAGILRQSFMYGAGYQVEKVISHPNYDSKTKNNDIALMKLQKPLTFNDLVKP VCLPNPGMMLQPEQLCWISGWGATEEKGKTSEVLNAAKVLLIETQRCNSRYVYD NLITPAMICAGFLQGNVDSCQGDSSGGLVTSKNNIWWLIGDTSWGS GCAKAYR PGVYGNVMVFTDWIYRQMRADG
Purity	Greater than 85% as determined by SDS-PAGE.
Source	E.coli
MW	53.7 kDa
Uniprot ID	O15393
Expiration Date	6 months from date of receipt.