

Product Datasheet

Rabbit Anti-Sheep IgG (H&L) antibody (AMCA) (orb705661)

Description	Rabbit polyclonal Secondary antibody against Rabbit Anti-Sheep IgG (H&L) conjugated to AMCA
Species/Host	Rabbit
Reactivity	Sheep
Conjugation	AMCA
Tested Applications	ELISA, FC, IF, IH
Immunogen	Sheep IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Dilution Range	IHC-P: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000
Expiration Date	12 months from date of receipt.