

## Product Datasheet

### Rabbit Anti-Sheep IgG (H&L) - AMCA (orb705661)

<b>Description</b>	Rabbit polyclonal Secondary antibody against Rabbit Anti-Sheep IgG (H&L) conjugated to AMCA
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Sheep
<b>Conjugation</b>	AMCA
<b>Tested Applications</b>	ELISA, FC, IF, IH
<b>Immunogen</b>	Sheep IgG
<b>Preservatives</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
<b>Form/Appearance</b>	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Rabbit
<b>Dilution Range</b>	IHC-P: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000
<b>Expiration Date</b>	12 months from date of receipt.

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)