

Product Datasheet

Goat Anti-Human IgG+IgM+IgA antibody (TRITC) (orb705680)

Description	Goat polyclonal Secondary antibody against Goat Anti-Human IgG+IgM+IgA conjugated to TRITC
Species/Host	Goat
Reactivity	Human
Conjugation	TRITC
Tested Applications	ELISA, FC, IF, IH
Immunogen	Human IgG+IgM+IgA
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	IHC-P: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000
Expiration Date	12 months from date of receipt.