

Product Datasheet

Rabbit Anti-Goat IgG (H&L) antibody (TRITC) (orb705692)

Description	Rabbit polyclonal Secondary antibody against Rabbit Anti-Goat IgG (H&L) conjugated to TRITC
Species/Host	Rabbit
Reactivity	Goat
Conjugation	TRITC
Tested Applications	ELISA, FC, IF, IH
Immunogen	Goat IgG
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Dilution Range	IHC-P: 1:100-1:1000, IF: 1:100-1:1000, FC: 1:1000-1:4000
Expiration Date	12 months from date of receipt.