

## Product Datasheet

### Myc-tag antibody (HRP) (orb705937)

<b>Description</b>	Rabbit polyclonal antibody against Anti-Myc-tag conjugated to HRP
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	WB
<b>Immunogen</b>	KLH-conjugated synthetic peptide encompassing a sequence of Myc-tag.
<b>Preservatives</b>	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Proclin
<b>Form/Appearance</b>	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Proclin
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Monoclonal
<b>Source</b>	Mouse
<b>Dilution Range</b>	WB: 1:2000-1:5000, ICC: 1:500-1:2000
<b>Expiration Date</b>	12 months from date of receipt.