

Product Datasheet

FLAG-tag antibody (Biotin) (orb705944)

Description	Rabbit polyclonal antibody against Anti-FLAG-tag conjugated to Biotin
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of FLAG-tag.
Preservatives	Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Source	Mouse
Dilution Range	WB: 1:2000-1:5000, ICC: 1:500-1:2000
Expiration Date	12 months from date of receipt.