

Product Datasheet

Beta-catenin (phospho-Y670) Blocking Peptide (orb705995)

Description	The peptide is used to block Anti-Beta-catenin (pY670) antibody < a href="https://www.biorbyt.com/ct..."
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CTNNB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9WU82 , Q02248 , P35222
Entrez	12387 , 84353 , 1499
Expiration Date	12 months from date of receipt.