

Product Datasheet

Caveolin 2 (phospho-Y19) Blocking Peptide (orb706039)

Description	The peptide is used to block Anti-Caveolin 2 (pY19) antibody < a href="https://www.biorbyt.com/cav2-...>
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CAV2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P51636 , Q9WVC3 , Q2IBC5
Entrez	363425 , 12390 , 858
Expiration Date	12 months from date of receipt.