

Product Datasheet

Acinus (phospho-S1180) Blocking Peptide (orb706046)

Description	The peptide is used to block Anti-Acinus (pS1180) antibody < a href="https://www.biorbyt.com/acin1-p..."
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ACIN1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9JIX8 , O95104 , Q63627 , A2AJT4 , Q9UKV3 , Q8VH51 , Q14498 , Q8TF01
Entrez	170791 , 57466 , 25957 , 22985 , 56215 , 9584
Expiration Date	12 months from date of receipt.