

Product Datasheet

TORC2 (phospho-S171) Blocking Peptide (orb706049)

Description	The peptide is used to block Anti-TORC2 (pS171) antibody < a href="https://www.biorbyt.com/crtc2-pho..."
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CRTC2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q3U182 , Q53ET0 , Q3LRZ1
Entrez	200186 , 74343 , 310615
Expiration Date	12 months from date of receipt.