

## Product Datasheet

# PKC beta (phospho-T641) Blocking Peptide (orb706051)

**Description** The peptide is used to block Anti-PKC beta (pT641) antibody < a href="https://www.biorbyt.com/prkcb-...

**Reactivity** Human, Mouse, Rat, Zebrafish

**Tested Applications** Blocking

**Target** PRKCB

**Form/Appearance** Lyophilized powder

**Storage** Shipped at 4°C. Store at -20°C for one year.

**Note** For research use only

**Application notes** Quality Control: The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.

**Purity** > 85%

**Source** Synthetic

**Uniprot ID** **P68404, P05771**

**Entrez** **5579, 18751**

**Expiration Date** 12 months from date of receipt.